

**Product Name :**  
Liquid Limit Device Hand Operated with Counter

**Product Code :**  
TEST5000RFQLAB-0001



**Description :**

Liquid Limit Device Hand Operated with Counter for Testing Lab Equipment Tender Suppliers India

**Technical Specification :**

The apparatus consists of a brass cup, a crank and cam mechanism, mounted on a rubber base of a specified hardness. A brass pin having a knurled end for easy removal holds the brass cup. The height of fall of the cup can be adjusted by a horizontal lead screw. The device is fitted with a drop counter to automatically record number of drops of the cup.