


**Product Name :**  
Bench type Consolidometer Front Loading Single Unit

**Product Code :**  
TEST5010RFQLAB-0001



**Description :**

Bench type Consolidometer Front Loading Single Unit for Testing Lab Equipment Tender Suppliers India

**Technical Specification :**

Loading Unit 20 kg/cm<sup>2</sup> capacity -1 No.  
Consolidation cell 60mm x 20mm height -1 No.  
Dial gauge 5 x 0.002mm -1 No.  
Set of Weights (Stress – up to 10kg/cm<sup>2</sup>) -1 Set