

**Product Name :**  
Drop Point of Lubricating Grease Apparatus Electrically  
Heated

**Product Code :**  
TEST5079RFQLAB-0004



**Description :**

Drop Point of Lubricating Grease Apparatus Electrically Heated

**Technical Specification :**

For determining drop point of lubricating grease.

It consists of a chromium plated brass cup 9.9 mm internal diameter at wider end, 40 ml beaker, metal rod, electric heater and stirrer. Suitable for operation on 220 V, 50 Hz single phase supply.