

**Product Name :**  
Laboratory CBR Apparatus Hand Operated

**Product Code :**  
TEST5011RFQLAB-0001



**Description :**

Laboratory CBR Apparatus Hand Operated for Testing Lab Equipment Tender Suppliers India

**Technical Specification :**

Hand Operated Load Frame, capacity 5000 kg (50 KN).  
CBR Mould 150 mm base plate and collar inner dia x 175 mm high.  
Cutting collar.  
Perforated swell plate.  
Spacer disc.  
Metallic ring made of mild steel.  
Metal tripod aluminium.  
Penetration Piston Assembly.  
Proving ring 10 kN capacity.  
One surcharge weight, 2.5 kg, annular.  
One surcharge weight, 2.5 kg slotted.  
Dial gauge 0.01 x 25mm

