

Product Name :
Bitumen Penetration Kit

Product Code :
TEST5080RFQLAB-0006



Description :

Bitumen Penetration Kit

Technical Specification :

Bitumen Penetration Kit for determining the consistency of semi-solid bituminous materials. It consists of a standard penetration needle mounted in a brass shank, a transfer dish 90 mm dia x 70 mm deep approx and two aluminium sample containers 55 mm dia x 35 mm deep and 55 mm dia x 57 mm deep.

We are a leading manufacturers of Bitumen Penetration Kit testing equipment for a range of industries, including materials science testing lab. We offer full range of Bitumen Penetration Kit testing equipment, including Bitumen Penetration Kit testing kit, Bitumen Penetration Kit testing machine, Bitumen Penetration Kit testing tools, Bitumen Penetration Kit test equipment, Bitumen Penetration Kit test machines, Bitumen Penetration Kit for Lab Testing Instruments. We are suppliers of materials testing tools and equipment for the construction industry. Our mission is to provide reliable, accurate, and easy-to-use Bitumen Penetration Kit testing equipment that helps customers build better, safer, and more sustainable infrastructure.

