

Product Name :
Transistorised Penetrometer Controller or Timer

Product Code :
TEST5080RFQLAB-0007



Description :

Transistorised Penetrometer Controller or Timer

Technical Specification :

This unit supplies 220V, D.C. output exactly for 5 seconds to a timing head (Built in Electromagnet). The timer head allows the plunger rod to fall freely for 5 seconds and clamps it automatically after 5 seconds. Suitable for operation on 220V, 50 Hz, single phase supply.

We are a leading manufacturers of Transistorised Penetrometer Controller or Timer testing equipment for a range of industries, including materials science testing lab. We offer full range of Transistorised Penetrometer Controller or Timer testing equipment, including Transistorised Penetrometer Controller or Timer testing kit, Transistorised Penetrometer Controller or Timer testing machine, Transistorised Penetrometer Controller or Timer testing tools, Transistorised Penetrometer Controller or Timer test equipment, Transistorised Penetrometer Controller or Timer test machines, Transistorised Penetrometer Controller or Timer for Lab Testing Instruments. We are suppliers of materials testing tools and equipment for the construction industry. Our mission is to provide reliable, accurate, and easy-to-use Transistorised Penetrometer Controller or Timer testing equipment that helps customers build better, safer, and more sustainable infrastructure.

Testing lab instruments India

