

Product Name :
Stop Watch Holder Bench Type

Product Code :
TEST5080RFQLAB-0010



Description :

Stop Watch Holder Bench Type

Technical Specification :

It is made of cast aluminium. The push button operates the watch for start, stop and reset by a single lever mechanism which is adjustable for different watch sizes. Supplied without stop watch.

We are a leading manufacturers of Stop Watch Holder Bench Type testing equipment for a range of industries, including materials science testing lab. We offer full range of Stop Watch Holder Bench Type testing equipment, including Stop Watch Holder Bench Type testing kit, Stop Watch Holder Bench Type testing machine, Stop Watch Holder Bench Type testing tools, Stop Watch Holder Bench Type test equipment, Stop Watch Holder Bench Type test machines, Stop Watch Holder Bench Type for Lab Testing Instruments. We are suppliers of materials testing tools and equipment for the construction industry. Our mission is to provide reliable, accurate, and easy-to-use Stop Watch Holder Bench Type testing equipment that helps customers build better, safer, and more sustainable infrastructure.

