

Product Name :
Unconfined Compression Apparatus Hand Operated
Proving Ring Type

Product Code :
TEST5012RFQLAB-0001



Description :

Unconfined Compression Apparatus Hand Operated Proving Ring Type for Testing Lab Equipment Tender Suppliers India

Technical Specification :

Screw jack with frame and dial gauge holder
Set of upper and lower platens, 150 mm dia.
Cone Seating – 2 Nos.
Proving Ring Adapter.
Proving ring fitted with dial gauge, 5 kN capacity.
Dial gauge 0.01 x 25 mm.