

Product Name :
Standard Tar Viscometer Electrically Heated

Product Code :
TEST5081RFQLAB-0002



Description :

Standard Tar Viscometer Electrically Heated for Testing Lab Equipment Tender Suppliers India

Technical Specification :

but the gas heating element is replaced by an immersion heating element and an auto transformer for controlling the temperature is provided.

Suitable for operation on 220V, 50 Hz, single phase supply.

Cup with 4 mm Orifice.

Viscometer Thermometers

| Type | Range |
|--------|--------------|
| IP 8C | 0° -44°C |
| IP 9C | 37.8° – 82°C |
| IP 10C | 76° – 122°C. |