

**Product Name :**  
Engler Viscometer Gas Heated

**Product Code :**  
TEST5081RFQLAB-0003



**Description :**

Engler Viscometer Gas Heated

**Technical Specification :**

For determining the viscosity of lubricating and fuel oils.

It has an oil cup with a lid carrying thermometer socket. The cup is fitted with a stainless steel jet, mounted in a highly polished spun brass water bath. The bath is fitted with a thermometer clip. The whole system is placed on the tripod stand with adjustable feet. Complete with a plug valve for the jet and a stirring device.

We are a leading manufacturers of Engler Viscometer Gas Heated testing equipment for a range of industries, including materials science testing lab. We offer full range of Engler Viscometer Gas Heated testing equipment, including Engler Viscometer Gas Heated testing kit, Engler Viscometer Gas Heated testing machine, Engler Viscometer Gas Heated testing tools, Engler Viscometer Gas Heated test equipment, Engler Viscometer Gas Heated test machines, Engler Viscometer Gas Heated for Lab Testing Instruments. We are suppliers of materials testing tools and equipment for the construction industry. Our mission is to provide reliable, accurate, and easy-to-use Engler Viscometer Gas Heated testing equipment that helps customers build better, safer, and more sustainable infrastructure.

---

# Testing lab instruments India

