

**Product Name :**  
Laboratory Vane Shear Apparatus Hand Operated

**Product Code :**  
TEST5013RFQLAB-0001



**Description :**

Laboratory Vane Shear Apparatus Hand Operated for Testing Lab Equipment Tender Suppliers India

**Technical Specification :**

A torque applicator having a base with a hole for holding the specimen mould in position, supporting frame carrying torque head adjustable in height and a graduated drum to measure the deformation of springs.

A set of four springs, one each of 2 kg cm, 4 kg cm, 6 kg cm and 8 kg cm.

A vane of size 12.0 x 24.0 mm long with a vane rod.

Specimen container to take a specimen of 50 mm dia x 75 mm high.

Wooden carrying case for the complete apparatus.