

**Product Name :**  
Automatic Blain Air Permeability Apparatus

**Product Code :**  
TEST5053RFQLAB-0003



**Description :**

Automatic Blain Air Permeability Apparatus

**Technical Specification :**

Automatic Blaine Apparatus is an automated version of Blaine apparatus conforming to testing procedure laid down in relevant National and International standards.

**Salient Features**

Fully PC controlled operation for testing and analysis

Automatic control of fluid movement

Automatic measurement of temperature & Viscosity during the test

Facility to configure and monitor various cement types

Easy and tamper-proof calibration using a lockable ball valve

Note- A PC is required for use with the Automatic Blaine Air Permeability Apparatus but it is not part of supply and has to be arranged by customer for loading software

**Supplied Complete with:**

Filter papers 12.8 mm dia. (1000 ea.)

Fill oil, 50 ml

Light grase, 50 ml

Plug, sensor internally, brush and funnel

With software, without PC

Power: 230V 50Hz, 110V/60Hz

We are a leading manufacturers of Automatic Blain Air Permeability Apparatus testing equipment for a range of

---

industries, including materials science testing lab. We offer full range of Automatic Blain Air Permeability Apparatus testing equipment, including Automatic Blain Air Permeability Apparatus testing kit, Automatic Blain Air Permeability Apparatus testing machine, Automatic Blain Air Permeability Apparatus testing tools, Automatic Blain Air Permeability Apparatus test equipment, Automatic Blain Air Permeability Apparatus test machines, Automatic Blain Air Permeability Apparatus for Lab Testing Instruments. We are suppliers of materials testing tools and equipment for the construction industry. Our mission is to provide reliable, accurate, and easy-to-use Automatic Blain Air Permeability Apparatus testing equipment that helps customers build better, safer, and more sustainable infrastructure.



## Testing lab instruments India