



Email : sales@testinglabinstruments.com

Product Name :
Shrinkage Bar Mould Single

Product Code :
TEST5056RFQLAB-0003



Description :

Shrinkage Bar Mould Single

Technical Specification :

For preparing specimens of 25 mm x 25 mm x 250 mm effective gauge length with inserted stainless steel reference studs for determining the soundness of Portland Cement by the autoclave test and potential alkali reactivity of cement-aggregate combination by mortar bar method. The mould consists of two side plates and two end plates assembled on a common base plate with screws and dowel pins. Two reference stud holders snugly fit into the mould and are held by screws through the end plates. Two knurled and threaded reference studs are fitted into these stud holders.

We are a leading manufacturers of Shrinkage Bar Mould Single testing equipment for a range of industries, including materials science testing lab. We offer full range of Shrinkage Bar Mould Single testing equipment, including Shrinkage Bar Mould Single testing kit, Shrinkage Bar Mould Single testing machine, Shrinkage Bar Mould Single testing tools, Shrinkage Bar Mould Single test equipment, Shrinkage Bar Mould Single test machines, Shrinkage Bar Mould Single for Lab Testing Instruments. We are suppliers of materials testing tools and equipment for the construction industry. Our mission is to provide reliable, accurate, and easy-to-use Shrinkage Bar Mould Single testing equipment that helps customers build better, safer, and more sustainable infrastructure.

Testing lab instruments India

