

Product Name :
Length Comparator

Product Code :
TEST5056RFQLAB-0007



Description :

Length Comparator

Technical Specification :

Used for measuring the length change of hardened cement paste, concrete and mortar and also due to potential alkali reactivity of cement aggregate combination in course of time.

The apparatus comprises the following :-

One frame consisting of one base plate, two pillars and a cross-bar. The cross bar is adjustable up and down on the threaded pillars.

One dial gauge, least count 0.002 mm, travel 5 mm which can be fixed to the cross bar. The plunger of the dial gauge has a concave spherical seating which sits over the reference studs of specimen bars while taking observation.

One reference bar of stainless steel with insulated handle.

We are a leading manufacturers of Length Comparator testing equipment for a range of industries, including materials science testing lab. We offer full range of Length Comparator testing equipment, including Length Comparator testing kit, Length Comparator testing machine, Length Comparator testing tools, Length Comparator test equipment, Length Comparator test machines, Length Comparator for Lab Testing Instruments. We are suppliers of materials testing tools and equipment for the construction industry. Our mission is to provide reliable, accurate, and easy-to-use Length Comparator testing equipment that helps customers build better, safer, and more sustainable infrastructure.



Testing lab instruments India