

**Product Name :**  
Triaxial Cell for 38/50mm diameter x 76/100mm High  
Specimen

**Product Code :**  
TEST5016RFQLAB-0001



**Description :**

Triaxial Cell for 38/50mm diameter x 76/100mm High Specimen for Testing Lab Equipment Tender Suppliers India

**Technical Specification :**

Suitable for lateral pressure up to 1000 kpa, the cell has four take off positions drilled in the base for cell pressure, pore water pressure, and for top and bottom drainage. Fitted with 4 no-volume change valves and one dead end plug. Cell is supplied complete with one loading pad, two porous stones, one sheath stretcher, one split sand former, 10 latex rubber sheaths and four sealing rings for 38/50 mm dia samples.

Except that the acrylic tube is reinforced for the cell to withstand a maximum lateral pressure of 2000kpa. It is complete with all accessories

