

**Product Name :**  
Briquette Mould Single

**Product Code :**  
TEST5057RFQLAB-0003



**Description :**

Briquette Mould Single

**Technical Specification :**

For casting of cement briquettes for tensile strength tests. It is a two part split mould made of gun metal. Two thumb screws facilitate easy and quick assembling and dismantling of the mould. The minimum cross section of the briquettes cast is 25.4 mm x 25.4 mm. Supplied complete with a steel base plate.

We are a leading manufacturers of Briquette Mould Single testing equipment for a range of industries, including materials science testing lab. We offer full range of Briquette Mould Single testing equipment, including Briquette Mould Single testing kit, Briquette Mould Single testing machine, Briquette Mould Single testing tools, Briquette Mould Single test equipment, Briquette Mould Single test machines, Briquette Mould Single for Lab Testing Instruments. We are suppliers of materials testing tools and equipment for the construction industry. Our mission is to provide reliable, accurate, and easy-to-use Briquette Mould Single testing equipment that helps customers build better, safer, and more sustainable infrastructure.

