

**Product Name :**  
Shrinkage Limit Apparatus

**Product Code :**  
TEST5001RFQLAB-0001



**Description :**

Shrinkage Limit Apparatus for Testing Lab Equipment Tender Suppliers India

**Technical Specification :**

75mm square Prong Plate made of 3mm thick acrylic sheet with three metal prongs.  
75mm square Plain Plate made of 3mm thick acrylic sheet.  
Stainless Steel Shrinkage Dish 45mm in dia, 15mm deep.  
Glass cup with ground edge, 50mm to 55mm dia, 25mm deep.  
Flexible Spatula with 80mm long, 20mm wide blade.

We are a leading manufacturers of Shrinkage Limit Apparatus testing equipment for a range of industries, including materials science testing lab. We offer full range of Shrinkage Limit Apparatus testing equipment, including Shrinkage Limit Apparatus testing kit, Shrinkage Limit Apparatus testing machine, Shrinkage Limit Apparatus testing tools, Shrinkage Limit Apparatus test equipment, Shrinkage Limit Apparatus test machines, Shrinkage Limit Apparatus testing Instruments. We are suppliers of materials testing tools equipment for the construction industry. Our mission is to provide reliable, accurate, and easy-to-use Shrinkage Limit Apparatus testing equipment that helps customers build better, safer, and more sustainable infrastructure.

---

# Testing lab instruments India

