

**Product Name :**  
Reaction Container

**Product Code :**  
TEST5067RFQLAB-0016



**Description :**

Reaction Container

**Technical Specification :**

The reaction container is one of the important apparatus, required for determination of the Potential Alkali Reactivity of aggregates by chemical method when used in concrete having high alkali cement. It is made from stainless steel and fitted with air tight cover, held in position with special clamp.

We are a leading manufacturers of Reaction Container testing equipment for a range of industries, including materials science testing lab. We offer full range of Reaction Container testing equipment, including Reaction Container testing kit, Reaction Container testing machine, Reaction Container testing tools, Reaction Container test equipment, Reaction Container test machines, Reaction Container for Lab Testing Instruments. We are suppliers of materials testing tools and equipment for the construction industry. Our mission is to provide reliable, accurate, and easy-to-use Reaction Container testing equipment that helps customers build better, safer, and more sustainable infrastructure.

