

Product Name :
Rock Cutting and Polishing Machine

Product Code :
TEST5034RFQLAB-0003



Description :

Rock Cutting and Polishing Machine for Testing Lab Equipment Tender Suppliers India

Technical Specification :

A single machine can perform both the operations i.e. cutting and polishing. The cutting size is limited upto NX core or smaller.

The machine is fixed on a pedestal. A sample holder is provided in which the sample is moved towards the cutting blade. After cutting the faces can be ground or polished to a very fine finish on the diamond impregnated polishing blade.

Technical specifications are as under :-

Dia of blade: 250/300 mm

Dia. of polishing disc : 200 mm

Operating HP: 2

Operated on: 440 Volts, 3 phase (220 volts Optional)

Water Flow valve : Provided

RPM of blade: 2200 (approx.)

Sample size Nx and height up to 150mm.

Jig for 150mm

Jig for 450mm

We are a leading manufacturers of Rock Cutting and Polishing Machine testing equipment for a range of

industries, including materials science testing lab. We offer full range of Rock Cutting and Polishing Machine testing equipment, including Cutting and Polishing Machine testing kit, Cutting and Polishing Machine testing machine, Rock Cutting and Polishing Machine testing tools, Cutting and Polishing Machine test equipment, Cutting and Polishing Machine test machines, Cutting and Polishing Machine for Lab Testing Instruments. We are suppliers of materials testing tools and equipment for the construction industry. Our mission is to provide reliable, accurate, and easy-to-use Rock Cutting and Polishing Machine testing equipment that helps customers build better, safer, and more sustainable infrastructure.



Testing lab instruments India