

**Product Name :**  
Hydraulic Power Supply

**Product Code :**  
TEST5045RFQLAB-0003



**Description :**

Hydraulic Power Supply

**Technical Specification :**

Hydraulic power supplies are compact in design and are suitable for the supply of required flow and pressure for the movement of the actuator. It has an oil tank of adequate capacity, vane type pump powered by a three phase motor. All the electrical controls including the temperature controller are fixed on one side of the tank. It includes all the accessories like pressure line filter, return line filter, oil level, relief valve, pressure gauge and shell & tube type heat exchanger. Anti vibration mountings are provided as standard along with the HPS. Cooling system for keeping oil within working temperature range is provided as standard. The system is kept at a distance from the loading unit and connected through flexible pipes.

Flow of the pump : 40LPM

Motor capacity : 25HP

Capacity of the tank : 200litres

Operating pressure : 200Bars

We are a leading manufacturer of Hydraulic Power Supply testing equipment for a range of industries, including materials science testing lab. We offer full range of Hydraulic Power Supply testing equipment, including Hydraulic Power Supply testing kit, Hydraulic Power Supply testing machine, Hydraulic Power Supply testing tools, Hydraulic Power Supply test equipment, Hydraulic Power Supply test machines, Hydraulic Power Supply for Lab Testing Instruments. We are suppliers of materials testing tools and equipment for the construction industry. Our mission is to provide reliable, accurate, and easy-to-use Hydraulic Power Supply testing equipment

---

that helps customers build better, safer, and more sustainable infrastructure.



## Testing lab instruments India