

Product Name :
Universal Permeameter

Product Code :
TEST5005RFQLAB-0002



Description :

Universal Permeameter for Testing Lab Equipment Tender Suppliers India

Technical Specification :

A wooden stand fitted with nine glass stand pipes of diameter ranging from 6 mm to 75 mm. A 75 mm dia tube is fitted with an overflow pipe. A rubber bung with brass tube is provided for constant head tests. Permeability mould 1000 ml capacity with detachable collar, 50 mm high, one base plate, metallic perforated swell plate and a solid false bottom plate for compaction of soil.

We are a leading manufacturers of Universal Permeameter testing equipment for a range of industries, including materials science testing lab. We offer full range of Universal Permeameter testing equipment, including Universal Permeameter testing kit, Universal Permeameter testing machine, Universal Permeameter testing tools, Universal Permeameter test equipment, Universal Permeameter test machines, Universal Permeameter testing Instruments. We are suppliers of materials testing tools equipment for the construction industry. Our mission is to provide reliable, accurate, and easy-to-use Universal Permeameter testing equipment that helps customers build better, safer, and more sustainable infrastructure.

