

**Product Name :**  
Flexible Wall Permeability Apparatus

**Product Code :**  
TEST5005RFQLAB-0005



**Description :**

Flexible Wall Permeability Apparatus

**Technical Specification :**

Permeability System As per ASTM D5084 The apparatus consists of a control panel with three lines each for In-fluent line Effluent line and Cell pressure line. Push in valves are provided in each line to set and measure pressure either on digital display or on master gauge. It has provision for additional modules with three lines for each cell. One module is for one cell.

The apparatus consists of a control panel with three lines each for

Master Control Panel 1

Additional Panel 1

Test chambers 1

De-airing system 1

We are a leading manufacturers of Flexible Wall Permeability Apparatus testing equipment for a range of industries, including materials science testing lab. We offer full range of Flexible Wall Permeability Apparatus testing equipment, including Flexible Wall Permeability Apparatus testing kit, Flexible Wall Permeability Apparatus testing machine, Flexible Wall Permeability Apparatus testing tools, Flexible Wall Permeability Apparatus test equipment, Flexible Wall Permeability Apparatus test machines, Flexible Wall Permeability Apparatus testing Instruments. We are suppliers of materials testing tools equipment for the construction industry. Our mission is to provide reliable, accurate, and easy-to-use Flexible Wall Permeability Apparatus testing equipment that helps customers build better, safer, and more sustainable infrastructure.



# Testing lab instruments India