



Email : sales@testinglabinstruments.com

**Product Name :**  
Pycnometer

**Product Code :**  
TEST5006RFQLAB-0005



**Description :**

Pycnometer

**Technical Specification :**

Test is performed as Suitable for determination of specific gravity of sand and fine gravel.  
The apparatus consists of a glass jar of 1 kg capacity, complete with brass cone, rubber seal and screwed cap.

We are a leading manufacturers of Pycnometer testing equipment for a range of industries, including materials science testing lab. We offer full range of Pycnometer testing equipment, including Pycnometer testing kit, Pycnometer testing machine, Pycnometer testing tools, Pycnometer test equipment, Pycnometer test machines, Pycnometer testing Instruments. We are suppliers of materials testing tools equipment for the construction industry. Our mission is to provide reliable, accurate, and easy-to-use Pycnometer testing equipment that helps customers build better, safer, and more sustainable infrastructure.



**Testing lab instruments India**