

Product Name :
High Speed Stirrer

Product Code :
TEST5007RFQLAB-0005



Description :

High Speed Stirrer

Technical Specification :

The test is performed as The unit comprises of :-

A mechanical stirring device, fitted with a motor which rotates a vertical haft having a stirring pedal at a speed of 8000 to 10,000 R.P.M.

Dispersion cup made of brass.

Baffle for use with dispersion cup.

We are a leading manufacturers of High Speed Stirrer testing equipment for a range of industries, including materials science testing lab. We offer full range of High Speed Stirrer testing equipment, including High Speed Stirrer testing kit, High Speed Stirrer testing machine, High Speed Stirrer testing tools, High Speed Stirrer test machines, High Speed Stirrer testing Instruments. We are suppliers of materials testing tools equipment for the construction industry. Our mission is to provide reliable, accurate, and easy-to-use High Speed Stirrer testing equipment that helps customers build better, safer, and more sustainable infrastructure.

Testing lab instruments India

