

Product Name :
Constant Volume Mould 38 mm dia x 126 mm Long

Product Code :
TEST5009RFQLAB-0008



Description :

Constant Volume Mould 38 mm dia x 126 mm Long, Constant volume Mould 50 mm dia x 180 mm long, Constant Volume Mould 75 mm dia x 260 mm long, Constant volume mould 100 mm x 350 mm long

Technical Specification :

Suitable for 38 mm dia x 76 mm long specimen, the mould is supplied complete with two 38 mm dia x 25 mm long plugs, split collar and 126 mm long ejecting plunger.

Suitable for a 50 mm dia x 100 mm long specimen, the mould is supplied with two 50 mm dia x 40 mm long end plugs, split collar and a 180 mm long ejecting plunger.

Suitable for a 75 mm dia x 150 mm long specimen, the mould is supplied complete with two 75 mm dia x 55 mm long end plugs, split collar and a 260 mm long ejecting plunger.

Suitable for 100 mm dia x 200 mm long specimen, the mould is supplied with two 100 mm dia x 75 mm long end plugs, split collar and a 350 mm long ejecting plunger.

