

**Product Name :**  
Three Gang Bench Type Consolidometer Front Loading

**Product Code :**  
TEST5010RFQLAB-0002



**Description :**

Three Gang Bench Type Consolidometer Front Loading for Testing Lab Equipment Tender Suppliers India

**Technical Specification :**

It is that the unit have 3 loading machines instead of one. The unit comprises:

Loading Unit 20 kg/cm<sup>2</sup> capacity – 3 No.

Consolidation cell 60mm x 20mm height – 3 No.

Dial gauge 5 x 0.002mm – 3 No.

Set of Weights (Stress – up to 10kg/cm<sup>2</sup>) – 3 Sets