

Product Name :
Three Gang Bench Type Consolidometer Front Loading

Product Code :
TEST5010RFQLAB-0002



Description :

Three Gang Bench Type Consolidometer Front Loading for Testing Lab Equipment Tender Suppliers India

Technical Specification :

It is that the unit have 3 loading machines instead of one. The unit comprises:

Loading Unit 20 kg/cm² capacity – 3 No.

Consolidation cell 60mm x 20mm height – 3 No.

Dial gauge 5 x 0.002mm – 3 No.

Set of Weights (Stress – up to 10kg/cm²) – 3 Sets