



Email : sales@testinglabinstruments.com

Product Name :
Unconfined Compression Apparatus Motorized Proving Ring Type

Product Code :
TEST5012RFQLAB-0002



Description :

Unconfined Compression Apparatus Motorized Proving Ring Type for Testing Lab Equipment Tender Suppliers India

Technical Specification :

Except that the loading unit is motorised and the gear system provides three different rates of strain 1.25 mm, 1.5 mm and 2.5 mm/min. The unit operates on 220V single phase supply.



Testing lab instruments India