

**Product Name :**  
Unsaturated Triaxial Cell

**Product Code :**  
TEST5017RFQLAB-0002



**Description :**

Unsaturated Triaxial Cell for Testing Lab Equipment Tender Suppliers India

**Technical Specification :**

This triaxial cell and universal sand testing machine is suitable for both unsaturated and saturated soil testing for soil samples of diameter 50mm & 100mm with L/D ratio of 2. It is a double walled Triaxial cell where inner wall have equal pressure on the both sides of the cell. It is suitable for confining pressure up to 1000kPa. The cell is provided with facility to measure both pore water pressure and pore air pressure during the test.

It is supplied complete with one loading pad, high air entry (HAE) value ceramic porous stones (2Bars) sealed on aluminium rings, two porous stones, one sheath stretcher, one split mould samplers, 10 latex rubber sheaths and four sealing rings for 50 & 100 mm diameter samples and drainage tubing.F

We are a leading manufacturers of Unsaturated Triaxial Cell testing equipment for a range of industries, including materials science testing lab. We offer full range of Unsaturated Triaxial Cell testing equipment, including Unsaturated Triaxial Cell testing kit, Unsaturated Triaxial Cell testing machine, Unsaturated Triaxial Cell testing tools, Unsaturated Triaxial Cell test equipment, Unsaturated Triaxial Cell test machines, Unsaturated Triaxial Cell for Lab Testing Instruments. We are suppliers of materials testing tools and equipment for the construction industry. Our mission is to provide reliable, accurate, and easy-to-use Unsaturated Triaxial Cell testing equipment that helps customers build better, safer, and more sustainable infrastructure.



# Testing lab instruments India