

**Product Name :**  
Electronic Volume Change Unit

**Product Code :**  
TEST5019RFQLAB-0004



**Description :**

Electronic Volume Change Unit

**Technical Specification :**

The unit is highly sensitive and accurate for the measurement of volume change. It has a frictionless rolling bellowfram seal coupled to submersible type displacement sensor. Sensitivity of the system is 0.01ml. Change in volume is directly displayed on digital display system in engineering units.

We are a leading manufacturers of Electronic Volume Change Unit testing equipment for a range of industries, including materials science testing lab. We offer full range of Electronic Volume Change Unit testing equipment, including Electronic Volume Change Unit testing kit, Electronic Volume Change Unit testing machine, Electronic Volume Change Unit testing tools, Electronic Volume Change Unit test equipment, Electronic Volume Change Unit test machines, Electronic Volume Change Unit for Lab Testing Instruments. We are suppliers of materials testing tools and equipment for the construction industry. Our mission is to provide reliable, accurate, and easy-to-use Electronic Volume Change Unit testing equipment that helps customers build better, safer, and more sustainable infrastructure.

