

Product Name :
Triaxial Cell With Accessories

Product Code :
TEST5020RFQLAB-0005



Description :

Triaxial Cell With Accessories

Technical Specification :

Standard Triaxial Cell can accommodate sample size, ranging from 38mm to 100mm diameter with L/D ratio of 1:2. It is suitable for both Static and Dynamic tests (Compression & Extension). It has a submersible load cell connected to the transfer bar. Linear bearings ensure smooth movement of the transfer bar (plunger). Base of the cell has four no-volume change valves, one each for back pressure, pore pressure, confining pressure and top drainage. Sensors for pore pressure and back pressure are attached to the stainless steel base.

Tests Possible : Static & Dynamic (Compression/Extension or both)

Confining Pressure : Up to 1000 kPa

Specimen Size : Up to 100mm diameter & 200mm height

Submersible Load cell : ± 5 kN

: ± 1 kN, ± 2 kN or ± 10 kN (Optional)

We are a leading manufacturers of Triaxial Cell With Accessories testing equipment for a range of industries, including materials science testing lab. We offer full range of Triaxial Cell With Accessories testing equipment, including Triaxial Cell With Accessories testing kit, Triaxial Cell With Accessories testing machine, Triaxial Cell With Accessories testing tools, Triaxial Cell With Accessories test equipment, Triaxial Cell With Accessories test machines, Triaxial Cell With Accessories for Lab Testing Instruments. We are suppliers of materials testing tools and equipment for the construction industry. Our mission is to provide reliable, accurate, and easy-to-use Triaxial Cell With Accessories testing equipment that helps customers build better, safer, and more sustainable

infrastructure.



Testing lab instruments India