

Product Name :
Triaxial Cell

Product Code :
TEST5020RFQLAB-0010



Description :

Triaxial Cell

Technical Specification :

Complete with adapters and piezoelectric sensor. Accommodates samples size of 50 mm diameter (optional 75mm & 100mm) with L:d ratio as 2:1. Confining pressure 1 MPa.
Cell is complete with Rubber Sheaths (Pack of 10) Porous Stones with central bore of 25 mm, sheath stretcher and O'rings (pack of 6).

We are a leading manufacturers of Triaxial Cell testing equipment for a range of industries, including materials science testing lab. We offer full range of Triaxial Cell testing equipment, including Triaxial Cell testing kit, Triaxial Cell testing machine, Triaxial Cell testing tools, Triaxial Cell test equipment, Triaxial Cell test machines, Triaxial Cell for Lab Testing Instruments. We are suppliers of materials testing tools and equipment for the construction industry. Our mission is to provide reliable, accurate, and easy-to-use Triaxial Cell testing equipment that helps customers build better, safer, and more sustainable infrastructure.

