

**Product Name :**  
Sieve Shaker Motorised

**Product Code :**  
TEST5023RFQLAB-0002



**Description :**

Sieve Shaker Motorised for Testing Lab Equipment Tender Suppliers India

**Technical Specification :**

Having a sturdy cast iron body, the shaker has an inclined sieve table which can accommodate a maximum of seven sieves of 200 mm dia. An adjustable top clamping plate is provided to hold the sieves. The table has a gyratory motion in addition to up and down jolting action. The operating gear assembly is enclosed in an oil bath fitted with a transparent oil level indicating window.

Suitable for operating on 220 v, 50 Hz, single phase supply.

We are a leading manufacturers of Sieve Shaker Motorised testing equipment for a range of industries, including materials science testing lab. We offer full range of Sieve Shaker Motorised testing equipment, including Sieve Shaker Motorised testing kit, Sieve Shaker Motorised testing machine, Sieve Shaker Motorised testing tools, Sieve Shaker Motorised test equipment, Sieve Shaker Motorised test machines, Sieve Shaker Motorised for Lab Testing Instruments. We are suppliers of materials testing tools and equipment for the construction industry. Our mission is to provide reliable, accurate, and easy-to-use Sieve Shaker Motorised testing equipment that helps customers build better, safer, and more sustainable infrastructure.

---

# Testing lab instruments India

