

Product Name :
Rolling Thin Film Oven Voltage 220V

Product Code :
TEST5111RFQLAB-0001



Description :

Rolling Thin Film Oven Voltage 220V for Testing Lab Equipment Tender Suppliers India

Technical Specification :

Mains voltage: 220V, 50HZ, 2.5kW
Control accuracy: 163 ± 0.5
Revolving speed of sample turntable: 5.5r/min
Time for each test: 5.25h
Dimension of the workroom: 430x410x450mm
Quantity of samples for each test: 4pcs