

**Product Name :**  
Rolling Thin Film Oven Voltage 220V

**Product Code :**  
TEST5111RFQLAB-0001



**Description :**

Rolling Thin Film Oven Voltage 220V for Testing Lab Equipment Tender Suppliers India

**Technical Specification :**

Mains voltage: 220V, 50HZ, 2.5kW  
Control accuracy:  $163 \pm 0.5$   
Revolving speed of sample turntable: 5.5r/min  
Time for each test: 5.25h  
Dimension of the workroom: 430×410×450mm  
Quantity of samples for each test: 4pcs