

Product Name :
Constant Temperature Water Bath 10L

Product Code :
TEST5114RFQLAB-0001



Description :

Constant Temperature Water Bath 10L for Testing Lab Equipment Tender Suppliers India

Technical Specification :

Constant temperature water bath is a laboratory instrument designed for factories, construction companies, colleges, universities, and research institutes. The water bath allows for multifunctional circulations, providing constant temperature for high or low temperature test.

Power supply: AC (220±10%)V, 50Hz

Bath volume: 10L

Temperature controlling range: -10~95

Temperature controlling accuracy: ±0.1

Circulating water: 3.5 L/min

Ambient temperature: -10~+30

Relative humidity: 85%

Dimension: 530x400x430mm (LxWxH)

Maximum power consumption: 1100W

We are a leading manufacturers of Constant Temperature Water Bath 10L testing equipment for a range of industries, including materials science testing lab. We offer full range of Constant Temperature Water Bath 10L testing equipment, including Constant Temperature Water Bath 10L testing kit, Constant Temperature Water Bath 10L testing machine, Constant Temperature Water Bath 10L testing tools, Constant Temperature Water Bath 10L test equipment, Constant Temperature Water Bath 10L test machines, Constant Temperature Water

Bath 10L for Lab Testing Instruments. We are suppliers of materials testing tools and equipment for the construction industry. Our mission is to provide reliable, accurate, and easy-to-use Constant Temperature Water Bath 10L testing equipment that helps customers build better, safer, and more sustainable infrastructure.



Testing lab instruments India