

Product Name :
Bitumen Mixture Centrifugal Extractor

Product Code :
TEST5118RFQLAB-0001



Description :

Bitumen Mixture Centrifugal Extractor for Testing Lab Equipment Tender Suppliers India

Technical Specification :

Bitumen mixture centrifugal extractor is designed and Standard Test Methods of Bitumen and Bituminous Mixtures for Highway Engineering.

This bitumen extractor is used to determine the percentage of bitumen in bituminous mixtures.

Feature

Stable operation with less vibration & noise;

Equip with speed control system, from minimum to maximum speed, easy to operate;

Standard 3000g bowl capacity and 1500g as optional.

Power supply: 220VAC 50HZ/60HZ

Motor power: 750W

Rotating speed: 3600rpm max.

Max. capacity: 3000g (1500g optional)

Overall dimension: Ø350mm x 526mm (H)

Weight: Around 70KG

Filter paper: 20pcs in standard configure

We are a leading manufacturers of Bitumen Mixture Centrifugal Extractor testing equipment for a range of industries, including materials science testing lab. We offer full range of Bitumen Mixture Centrifugal Extractor testing equipment, including Bitumen Mixture Centrifugal Extractor testing kit, Bitumen Mixture Centrifugal

Extractor testing machine, Bitumen Mixture Centrifugal Extractor testing tools, Bitumen Mixture Centrifugal Extractor test equipment, Bitumen Mixture Centrifugal Extractor test machines, Bitumen Mixture Centrifugal Extractor for Lab Testing Instruments. We are suppliers of materials testing tools and equipment for the construction industry. Our mission is to provide reliable, accurate, and easy-to-use Bitumen Mixture Centrifugal Extractor testing equipment that helps customers build better, safer, and more sustainable infrastructure.



Testing lab instruments India