

Product Name :
Paraffin Content Tester Machine

Product Code :
TEST5120RFQLAB-0001



Description :

Paraffin Content Tester Machine for Testing Lab Equipment Tender Suppliers India

Technical Specification :

Paraffin content tester is suitable for determining the Paraffin content in the petroleum asphalt.

Power supply: AC (220±10%)V, 50Hz

Heating power: 700W

Cooling power: 1000W

Stirring rate of motor: 1200r/min

Temperature sensor: Pt100

Refrigeration medium: F22

Temperature accuracy: ±0.1°

Working environment: 30°~40°, RH?85%

Total power consumption: ?1800W

Overall dimension: 620mm×500mm×1600mm (The air flow distributor is included)

We are a leading manufacturers of Paraffin Content testing equipment for a range of industries, including materials science testing lab. We offer full range of Paraffin Content testing equipment, including Paraffin Content testing kit, Paraffin Content testing machine, Paraffin Content testing tools, Paraffin Content test equipment, Paraffin Content test machines, Paraffin Content for Lab Testing Instruments. We are suppliers of materials testing tools and equipment for the construction industry. Our mission is to provide reliable, accurate, and easy-to-use Paraffin Content testing equipment that helps customers build better, safer, and more sustainable infrastructure.



Testing lab instruments India