

Product Name :
Aggregate Impact Value Testing Apparatus Machine

Product Code :
TEST5128RFQLAB-0001



Description :

Aggregate Impact Value Testing Machine for Testing Lab Equipment Tender Suppliers India

Technical Specification :

Impact hammer weight: 13.75kg±0.05kg
Hammer falling height: 380±5mm
Impacting cup: ?102×50mm
Test sieve: 15.0mm, 10.0mm, 2.50mm
Metal tamping rod with one end in hemispherical shape: ?10×230mm
Metal measuring cylinder: ?75±1×50±1mm

We are a leading manufacturers of Aggregate Impact Value Testing Apparatus testing equipment for a range of industries, including materials science testing lab. We offer full range of Aggregate Impact Value Testing Apparatus testing equipment, including Aggregate Impact Value Testing Apparatus testing kit, Aggregate Impact Value Testing Apparatus testing machine, Aggregate Impact Value Testing Apparatus testing tools, Aggregate Impact Value Testing Apparatus test equipment, Aggregate Impact Value Testing Apparatus test machines, Aggregate Impact Value Testing Apparatus for Lab Testing Instruments. We are suppliers of materials testing tools and equipment for the construction industry. Our mission is to provide reliable, accurate, and easy-to-use Aggregate Impact Value Testing Apparatus testing equipment that helps customers build better, safer, and more sustainable infrastructure.



Testing lab instruments India