



Email : sales@testinglabinstruments.com

**Product Name :**  
Sand Equivalent Test Apparatus

**Product Code :**  
TEST5129RFQLAB-0001



**Description :**

Sand Equivalent Test Apparatus for Testing Lab Equipment Tender Suppliers India

**Technical Specification :**

- Rated power: 160W
- Vibrating frequency: 90 times/min
- Amplitude: 203mm±1.0mm
- Self-control time: 30s/1 time
- Accessory box
- Test tube: 4pcs, inner dia. (32±0.25mm), height (420±0.25mm)
- Piston: 1pc
- Washing pipe: 1pc
- Thick liquor: 8bottles

We are a leading manufacturers of Sand Equivalent Test Set testing equipment for a range of industries, including materials science testing lab. We offer full range of Sand Equivalent Test Set testing equipment, including Sand Equivalent Test Set testing kit, Sand Equivalent Test Set testing machine, Sand Equivalent Test Set testing tools, Sand Equivalent Test Set test equipment, Sand Equivalent Test Set test machines, Sand Equivalent Test Set for Lab Testing Instruments. We are suppliers of materials testing tools and equipment for the construction industry. Our mission is to provide reliable, accurate, and easy-to-use Sand Equivalent Test Set testing equipment that helps customers build better, safer, and more sustainable infrastructure.



# Testing lab instruments India