

Product Name :
Bitumen Automatic Softening Point Apparatus Ring and Ball

Product Code :
TEST5109RFQLAB-0002



Description :

Bitumen Automatic Softening Point Apparatus Ring and Ball for Testing Lab Equipment Tender Suppliers India

Technical Specification :

Automatic softening point apparatus is used to determine softening point of petroleum asphalt, coal pitch, liquid petroleum asphalt, and various asphalt products. It's widely used for asphalt manufacturing companies, construction companies, colleges & universities, scientific research institutes, etc.

Test medium: Water or glycerin (the instrument can automatically distinguish)

Temperature measurement range: -50.0?-125.0?

Test temperature range: 5?-125?

Temperature control rate: 5? \pm 0.5?/min

Test time: 00 minutes 00 seconds-99 min 59 seconds

Temperature control method: Automatic control, magnetic stirring

Steel ball weight:

Steel ball diameter: 3.50 \pm 0.05g

9.53mm

Beaker: 1000ml

Heater: 700W

Testing lab instruments India

