



Email : sales@testinglabinstruments.com

Product Name :
Bitumen Ductility Machine 4500W

Product Code :
TEST5110RFQLAB-0002



Description :

Bitumen Ductility Machine 4500W for Testing Lab Equipment Tender Suppliers India

Technical Specification :

Ductility testing machine is intended for determining the length at which the bitumen specimen is elongated before breaking at a specific speed and temperature.

Power: AC (220±10%)V, 50Hz

Measurement length: 1.5m (±10mm)

Heating power: 4500W

Temperature control: Range: 5~35?, accuracy: ±0.1?

Traction speed: 10~50mm/min, stepless adjustment, accuracy: ±1mm

Traction force: Range: 0~300N, accuracy: ±1N, resolution: 0.1N

Data output: (1) Communication port RS232 (2) Printed by micro-printer

Refrigeration: Compressor 1.0P, input power is 800W

Working temperature -10? ~+35?

Relative humidity: ?85%

Maximum power consumption: 5600W

Testing lab instruments India

