

Product Name :
Rolling Thin Film Oven Testing Equipment

Product Code :
TEST5111RFQLAB-0002



Description :

Rolling Thin Film Oven Testing Equipment for Testing Lab Equipment Tender Suppliers India

Technical Specification :

Working power: AC (220±10%)V, 50Hz
Power consumption: 2.6kW
Working temperature: 163°C
Temperature control accuracy: ±0.5°C
Plate rotation speed: (15±0.2)r/min
Air flow: (4000±200)mL/min
Timing device: alarm at 85 minutes
Working temperature: 5°C~50°C
Relative humidity: ?85%

We are a leading manufacturers of Rolling Thin Film Oven testing equipment for a range of industries, including materials science testing lab. We offer full range of Rolling Thin Film Oven testing equipment, including Rolling Thin Film Oven testing kit, Rolling Thin Film Oven testing machine, Rolling Thin Film Oven testing tools, Rolling Thin Film Oven test equipment, Rolling Thin Film Oven test machines, Rolling Thin Film Oven for Lab Testing Instruments. We are suppliers of materials testing tools and equipment for the construction industry. Our mission is to provide reliable, accurate, and easy-to-use Rolling Thin Film Oven testing equipment that helps customers build better, safer, and more sustainable infrastructure.



Testing lab instruments India