

**Product Name :**  
Automatic Mixture Blender

**Product Code :**  
TEST5115RFQLAB-0002



**Description :**

Automatic Mixture Blender for Testing Lab Equipment Tender Suppliers India

**Technical Specification :**

Automatic mixture blender adopts advanced microprocessor controlling technology to present the blending temperature ranging from ambient temperature to 200?. This automatic mixture blender can automatically control PID parameters and blending time, which is suitable to be used in the laboratory for the sample test.

Blending volume: 60L

Temperature control range: Room temperature~200?

Temperature control accuracy:  $\pm 5^\circ$

Timing range: 0~999s

Timing accuracy:  $\pm 0.1$ s

Rotation rate of paddle: Rotation speed: 67 r/min Revolution speed: 42 r/min

Blending motor: 380V, 50Hz, 2200W, 1400r/min

Lifting motor: 380V, 50Hz, 750W, 1400 r/min

Overall dimension: 1400mmx750mmx2400mm

Power supply: AC (380 $\pm$ 10%)V, 50Hz, 30A (three-phase four-wire system, with zero line)

Ambient temperature: -5?~+50?

Relative humidity:  $\leq 80\%$

We are a leading manufacturers of Automatic Mixture Blender testing equipment for a range of industries, including materials science testing lab. We offer full range of Automatic Mixture Blender testing equipment, including Automatic Mixture Blender testing kit, Automatic Mixture Blender testing machine, Automatic Mixture

---

Blender testing tools, Automatic Mixture Blender test equipment, Automatic Mixture Blender test machines, Automatic Mixture Blender for Lab Testing Instruments. We are suppliers of materials testing tools and equipment for the construction industry. Our mission is to provide reliable, accurate, and easy-to-use Automatic Mixture Blender testing equipment that helps customers build better, safer, and more sustainable infrastructure.



## Testing lab instruments India