

**Product Name :**  
Binder Extractor Apparatus

**Product Code :**  
TEST5118RFQLAB-0003



**Description :**

Binder Extractor Apparatus

**Technical Specification :**

Centrifugal extractor is used to perform reliable analysis on bituminous mixtures for quantitative determination of binder or bitumen contained in pavement samples and hot mixed mixtures.

Capacity: 1000g~1500g

Extracting precision:  $\pm 0.1\%$

Extracting time: (20~40)minutes per time

Centrifuging speed: 5500RPM, 11000RPM

Working environment: 5 $^{\circ}$ ~40 $^{\circ}$ , RH $\leq 80\%$

Power supply: AC (380V $\pm 10\%$ ), 50 Hz

Maximum power consumption: About 5kW

Pressure of cooling water:  $\leq 2$ bar

Temperature of cooling water:  $\leq 12^{\circ}$

We are a leading manufacturers of Binder Extractor testing equipment for a range of industries, including materials science testing lab. We offer full range of Binder Extractor testing equipment, including Binder Extractor testing kit, Binder Extractor testing machine, Binder Extractor testing tools, Binder Extractor test equipment, Binder Extractor test machines, Binder Extractor for Lab Testing Instruments. We are suppliers of materials testing tools and equipment for the construction industry. Our mission is to provide reliable, accurate, and easy-to-use Binder Extractor testing equipment that helps customers build better, safer, and more sustainable infrastructure.



# Testing lab instruments India