


**Product Name :**  
Los Angeles Abrasion Tester Machine

**Product Code :**  
TEST5136RFQLAB-0001



**Description :**

Los Angeles Abrasion Tester Machine for Testing Lab Equipment Tender Suppliers India

**Technical Specification :**

Revolution: 30–33r/min  
Total mass of steel balls: 5000±5g/12pcs  
Mass of each steel ball: 390–445g  
Power supply: AC 220V, 1.5kW

We are a leading manufacturers of Los Angeles Abrasion testing equipment for a range of industries, including materials science testing lab. We offer full range of Los Angeles Abrasion testing equipment, including Los Angeles Abrasion testing kit, Los Angeles Abrasion testing machine, Los Angeles Abrasion testing tools, Los Angeles Abrasion test equipment, Los Angeles Abrasion test machines, Los Angeles Abrasion for Lab Testing Instruments. We are suppliers of materials testing tools and equipment for the construction industry. Our mission is to provide reliable, accurate, and easy-to-use Los Angeles Abrasion testing equipment that helps customers build better, safer, and more sustainable infrastructure.

