

**Product Name :**  
Freeze Thaw Cabinet Antifreeze Fluid

**Product Code :**  
TEST5139RFQLAB-0002



**Description :**

Freeze Thaw Cabinet Antifreeze Fluid for Testing Lab Equipment Tender Suppliers India

**Technical Specification :**

Freeze thaw cabinet is suitable for the determination of concrete specimen under the water freezing and thawing conditions. This freezing and thawing testing machine measures the resistance of concrete to deterioration caused by repeated cycles of freezing and thawing in order to evaluate concrete performance.

Accepted specimen: 10pcs

Specimen dimension: 100×100×400mm

Freezing and thawing circle: Around 2.5-4 hours

While thawing time is no less than 1/4 of the whole circle

Freezing and thawing medium: Antifreeze fluid (-45?)

Cooling compressor: 3HP

Heating power: 2.5kW

Max. working power: 3kW

Power supply: 3N-380V, 50Hz

We are a leading manufacturers of Freeze Thaw Cabinet Antifreeze Fluid testing equipment for a range of industries, including materials science testing lab. We offer full range of Freeze Thaw Cabinet Antifreeze Fluid testing equipment, including Freeze Thaw Cabinet Antifreeze Fluid testing kit, Freeze Thaw Cabinet Antifreeze Fluid testing machine, Freeze Thaw Cabinet Antifreeze Fluid testing tools, Freeze Thaw Cabinet Antifreeze Fluid test equipment, Freeze Thaw Cabinet Antifreeze Fluid test machines, Freeze Thaw Cabinet Antifreeze Fluid for Lab Testing Instruments. We are suppliers of materials testing tools and equipment for the construction industry.

---

Our mission is to provide reliable, accurate, and easy-to-use Freeze Thaw Cabinet Antifreeze Fluid testing equipment that helps customers build better, safer, and more sustainable infrastructure.



## Testing lab instruments India