

**Product Name :**  
Flexure Testing Equipment

**Product Code :**  
TEST5144RFQLAB-0005



**Description :**

Flexure Testing Equipment

**Technical Specification :**

Flexure testing machine is intended for flexural strength test of cement specimen and other non-metallic material.

Max. loading capacity: Single lever: 1000N, double-lever: 5000N

Measuring accuracy:  $\pm 1\%$

Single lever output ratio: 10:1

Double-lever output ratio: 50:1

Loading speed ratio of double-lever: 50N/s

Support roller distance:  $100 \pm 0.1\text{mm}$

Baffle plate span: 46mm

We are a leading manufacturers of Flexure testing equipment for a range of industries, including materials science testing lab. We offer full range of Flexure testing equipment, including Flexure testing kit, Flexure testing machine, Flexure testing tools, Flexure test equipment, Flexure test machines, Flexure for Lab Testing Instruments. We are suppliers of materials testing tools and equipment for the construction industry. Our mission is to provide reliable, accurate, and easy-to-use Flexure testing equipment that helps customers build better, safer, and more sustainable infrastructure.



**Testing lab instruments India**