

**Product Name :**  
Soundness Steam tank

**Product Code :**  
TEST5056RFQLAB-0013



**Description :**

Soundness Steam tank

**Technical Specification :**

Used for the determination of the soundness of building limes, use with le-chatelier-moulds.

Technical Parameters

Max. temperature: 100?

Nominal volume: 31L

Time for water temperature rise (from 20? to 100?): 30±5min

Heating time: 0-3.5 h

Heater power: 220VAC 50Hz 4kW

We are a leading manufacturers of Soundness Steam tank testing equipment for a range of industries, including materials science testing lab. We offer full range of Soundness Steam tank testing equipment, including Soundness Steam tank testing kit, Soundness Steam tank testing machine, Soundness Steam tank testing tools, Soundness Steam tank test equipment, Soundness Steam tank test machines, Soundness Steam tank for Lab Testing Instruments. We are suppliers of materials testing tools and equipment for the construction industry. Our mission is to provide reliable, accurate, and easy-to-use Soundness Steam tank testing equipment that helps customers build better, safer, and more sustainable infrastructure.

---

# Testing lab instruments India

