

Product Name :
Cement Materials Test Mill

Product Code :
TEST5146RFQLAB-0001



Description :

Cement Materials Test Mill for Testing Lab Equipment Tender Suppliers India

Technical Specification :

Cement test mill is used to grind cement clinker and other materials in the laboratory.

Inside dia. of barrel and length: ?500x500mm

Grinding medium: 100kg

Material size: Loading capacity: ?5kg

Barrel speed: ?48r/min

Grinding time: ?30min

We are a leading manufacturers of Cement Materials Test Mill testing equipment for a range of industries, including materials science testing lab. We offer full range of Cement Materials Test Mill testing equipment, including Cement Materials Test Mill testing kit, Cement Materials Test Mill testing machine, Cement Materials Test Mill testing tools, Cement Materials Test Mill test equipment, Cement Materials Test Mill test machines, Cement Materials Test Mill for Lab Testing Instruments. We are suppliers of materials testing tools and equipment for the construction industry. Our mission is to provide reliable, accurate, and easy-to-use Cement Materials Test Mill testing equipment that helps customers build better, safer, and more sustainable infrastructure.

Testing lab instruments India

